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ELECTRODES – MANUAL / EN

ELEKTRODEN – ANLEITUNG / DE

ELECTRODES – INSTRUCTIONS / FR

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EN / User manual pH electrodes

01. Check electrode optically for damages.
02. Remove watering cap carefully by unscrewing the lid and flush electrode with water. Carefully dry electrode with a soft tissue. Do not rub. Rubbing can delay the response time.
03. Check interior of the electrode glass membrane at the tip for air bubbles. Remove existing air bubbles through shaking the electrode.
04. The electrode is now ready to use. Perform calibration and measurement according to the user manual of the pH-meter. Always use fresh buffer solution, the alkaline and neutral buffer can pick up CO₂ from the air, what can lead to a change in pH-value.
05. Dry stored electrodes show instable values. For regeneration, store the electrode overnight in 3M KCl-solution (23 g KCl in 100 ml distilled water). Always keep hP-electrodes in humid conditions.
06. To avoid loss of KCl from the gel-electrolyte through diffusion, the electrodes should only be stored in 3M KCl-solution.
07. The life expectancy of pH-electrodes can vary a lot through measurement conditions like temperature and composition of the measurement solution. In very mild conditions the life expectancy can be up to 3 years at room temperature and 1-3 months at 90°C. A slow ageing process also shows during storage time therefore it is limited.
08. pH electrodes with plastic shaft have limited chemical resistance. Avoid long term contact with concentrated acids and bases, alcohols, aromatic hydrocarbons, esters and ketones.

EN / User manual Redox electrode

For Redox (ORP-) electrodes basically the same rules apply like to the pH-electrodes. There are commercially obtainable redox-buffers with ca. 465 mV that serve as function check of the electrode. (The redox electrode cannot be calibrated like a pH electrode).

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